



E SEI IN MANI SICURE!

REFLEXX S.r.l.

INFORMATIVE NOTE

rev 01 del 04 maggio 2015

Product: reflexx 97 art. R97/ M - art. R97/ L - art. R97/XL - art. R97/ XXL

Emessa da G.Isetti - Presidente

Pag.: 1/4

NITRILE GLOVES FLOCKLINED GREEN TRADEMARK REFLEXX 97

1. NAME AND ADDRESS OF MANUFACTURER:

REFLEXX S.R.L.

Via Passeri, 2

46019 VIADANA (MN)

Tel. 0375 833164 Fax. 0375 464504

2. GENERAL DESCRIPTION OF PRODUCT:

- nitrile synthetic glove with five fingers, flocklined, green color, textured fingertips and antiskid surface
- non-sterile gloves, powder free , internally chlorinated
- multi-purpose glove to protect the hand and / or the product and can be used for general purposes.
- glove used as a medical device by examination, even invasive , intended for temporary use
- glove can be used as Personal Protective Equipment (PPE) also against chemicals and micro-organisms
- Unit Sale in pairs in envelope
- Outer packagings containing No. 50 envelopes

Imballo		
Cod./Taglie	Confezione	Cartone
R97/M	8032891632499	8032891637494
R97/L	8032891632505	8032891637500
R97/XL	8032891632512	8032891637517
R97/XXL	8032891632529	8032891637524

3. QUALITY AND COMPLIANCE WITH REGULATIONS:

- As a PPE category III risk according to IT Legislative Decree 475/92 implementing Dir 89/686/EEC (SATRA _0321). The product is certified according to the EN 420, EN 374 I, II e III e EN 388.
- The product may come into contact with food as it complies with the relevant legislation: IT Directives 82/711/EEC, 85/572/EEC, 93/8/EEC, 97/48/EC, 2002/72/EC and subsequent amendments, regulations 1935/2004/EC, 1895/2005/CE, EU n.10/2011, DM 21/03/1973 and subsequent updates and changes, Presidential Decree 777/82 and subsequent updates and changes in accordance with the relevant regulations UNI EN 1186 1-15:2003. For such use, follow the instructions for use.



E SEI IN MANI SICURE!

REFLEXX S.r.l.	INFORMATIVE NOTE	rev 01 del 04 maggio 2015
Product: reflexx 97 art. R97/ M - art. R97/ L - art. R97/XL - art. R97/ XXL		
<i>Emessa da G.Isetti - Presidente</i>	Pag.: 2/4	

4. DIMENSIONS AND PERFORMANCE:

Dimensions (in mm.)	size		Length (min.) (± 10)	Thickness (mm) (± 0.05)
	description size	size CE		
	M	8	330	0.40
	L	9	330	
	XL	10	330	
	XXL	11	330	

EN 420	construction compliant	in compliance
	measures	in compliance
	donning	in compliance

EN 374	test water loss	passed
	evidence of air leak	passed

EN 374-3	Performance level of permeation	
	methanol	2
	Sodium hydroxide 50%	6
	N-heptane	6

EN 388		Protection Levels
	Resistance to abrasion	4
	Resistance to blade cut and tear	0
	Resistance to perforation	2

EN 455-2:2009/A1	Tensile strength	≥6 N. (min)
------------------	------------------	-------------

For the certification of PPE category III, gloves reflexx 97 are tested with the following chemicals (EN 374): methanol,



E SEI IN MANI SICURE!

REFLEXX S.r.l.

INFORMATIVE NOTE

rev 01 del 04 maggio 2015

Product: reflexx 97 art. R97/ M - art. R97/ L - art. R97/XL - art. R97/ XXL

Emessa da G.Isetti - Presidente

Pag.: 3/4

Sodium hydroxide 50%, N-heptane and are tested (legislation EN-388) for resistance to abrasion, to cutting by blade, for tearing and perforation.

Gloves reflexx 97 are controlled according to the sampling method ISO 2859 AQL 1.5 (level G1) for freedom from holes.

5. INSTRUCTIONS FOR USE AND STORAGE

1. Please read carefully the information note before using the glove
2. Please keep this information note for the duration of the useful life of the glove
3. Disposable nitrile glove - PPE : Risk Category III can also be used against chemicals and micro-organisms (Directive 89/686/EEC)
4. Based on the tests carried out the glove Reflexx 97 is suitable for contact with aqueous, acid, alcoholic, dairy, oily or fatty, and dried foods for which it is necessary the use of simulants A, B, C, D1 and D2 for 30 minutes at 40°C.

The glove may come into contact with all types of foods, for example with non-alcoholic beverages and alcoholic, bakers and pastry wares fatty substances on the surface; Chocolate, sugar and products thereof Confectionery fatty substances on the surface; Whole fruit, fresh or chilled, unpeeled, nuts (eg chestnuts, peanuts, almonds and the like), also roasted, peeled and in the form of pasta and cream; Whole vegetables, fresh or chilled, unpeeled and fresh vegetables, peeled or in pieces; fats and animal and vegetable oils (margarine, butter and other fats and oils made from water emulsions in oil); fish, fresh, chilled, salted and smoked including fish eggs, canned fish; crustaceans and molluscs; meats of all kinds, fresh, chilled, salted, smoked, processed meat products (eg ham, salami) or in the form of paste or cream, produced meat marinated in an oily medium, canned meat; whole eggs, egg yolk, egg white, liquid and cooked; milk and milk-based drinks, yoghurts and similar natural cheeses and cheese fondants, soft cheeses and cheese preserved in aqueous or oily medium (mozzarella, feta, etc.); vinegar; foods fried or roasted; sauces (mayonnaise and other oil / water emulsions, fruit pickles); sandwich, toast, ice cream, frozen food and frozen food; Cocoa paste, spices and herbs in an oily medium (pesto, curry paste). For contact with specific foods and / or special conditions of use require the supplier information to the specific conditions of use.

5. The product does not contain natural rubber latex
6. In rare cases, you may encounter transient skin reactions and hypersensitivity reactions. In case of a reaction discontinue use . Nitrile gloves Reflexx 97 are designed in order to eliminate allergic reactions associated with latex proteins . The absence of dust enables the use of the glove even to people who suffer from irritations caused by allergies to corn starch or a particularly sensitive skin
7. Avoid contact with aggressive substances such as fats , oils , acids, hydrocarbons and other chemical agents ; has low resistance to ketone solvents.



E SEI IN MANI SICURE!

REFLEXX S.r.l.

INFORMATIVE NOTE

rev 01 del 04 maggio 2015

Product: reflexx 97 art. R97/ M - art. R97/ L - art. R97/XL - art. R97/ XXL

Emessa da G.Isetti - Presidente

Pag.: 4/4

8. Before using the product, test it in real conditions of use
9. The glove is a device subject to wear and tear so it is recommended the frequent replacement , and in all cases where visually you can see imperfections. Always wear gloves your hands are dry and clean
10. Keep the packaging for further information and to ensure traceability
11. Keep the product in its original packaging in a dry place at a temperature not exceeding 40 ° C or below -10 ° C
12. Avoid exposure to sunlight and heat sources
13. The expiry date is indicated on the packaging and it is valid when the product must is properly stored

REFLEXX Srl

Via Passeri, 2 - 46019 VIADANA (MN)

P.I. e C.F. 02085450209

Tel. 0375.833164 - Fax 0375.464504

e-mail: info@reflexx.com