


INFORMATIVE NOTE
Name: REFLEXX 70
Date: 15/05/2024 Rev: 021
MANUFACTURER

 REFLEXX SpA
 Via Passeri, 2 46019 VIADANA (MN)
 Tel.: + 39 0375 758891 Fax + 39 0375
 464504 info@reflexx.com
www.reflexx.com

Pack

Cod./Size	Box	Carton
R70/ S (7)	8032891631485	8032891636480
R70/ M (8)	8032891631492	8032891636497
R70/ L (9)	8032891631508	8032891636503
R70/ XL (10)	8032891631515	8032891636510
R70/ XXL (11)	8032891633519	8032891638514

DESCRIPTION: ■ disposable examination glove powder free, not sterile ■ **Material:** nitrile ■ **Shape:** ambidextrous; five fingers, beaded cuff ■ **Colour:** blue ■ **Feature:** internally chlorinated; textured surface and textured fingertips ■ **Packing:** 100 pieces/box; 10 boxes/carton ■ **ISO 2859:** AQL 1.5 (level G1) freedom from holes; AQL 2.5 (level G1) major defects; AQL 4 (level G1) minor defects

CONFORMITY AND PERFORMANCES
Personal Protective Equipment cat III

0465 ANCI_CIMAC, Via Aguzzafame, 60/b 27029 Vigevano (PV)

Reg EU 2016/425 - EN ISO 374-1:2016 TYPE B - EN ISO 374-5:2016 - EN 374-2:2014 - EN 420:2003 + A1 2009

Medical Device

Reg EU 2017/745 – EN 455 1,2,3 & 4

(CND T01020204 NON SURGICAL NITRILE GLOVES)



Reg CE 1935/2004 – Reg EU 10/2011 - DM 21/03/73. To such use carefully read the declaration of suitability for food contact (DICA).

DIMENSIONS



Cod./Size	Length (mm)	Width mm (±5)
R70/ S (7)	230*	85
R70/ M (8)	240	95
R70/ L (9)	250	105
R70/ XL (10)	260	115
R70/ XXL (11)	270	125

* Length min. 240 mm

Thickness mm	Palm (+/- 0.02)	0.10
------------------------	-----------------	------

EN 455-2:2009/A1	Tensile strength	≥6 N
-------------------------	------------------	------

EN ISO 374-5:2016	Protection from microorganisms (bacteria and fungi)	Virus protection ISO 16604:2004
	pass	pass


EN ISO 374-1:2016

Chemical compound:	COD	EN 16523-1:2015 TYPE B		EN 374-4:2013
		Performance level	Breakthrough time	% degradation
Sodium hydroxide 40%	K	3	>60'	41,3%
Hydrogen peroxide 30%	P	2	>30'	70,6%
Formaldehyde 37%	T	2	>30'	56,3%

RISKS FROM WHICH THE DPI IS DESIGNED TO PROTECT:

COD		Protection
K Sodium hydroxide 40%	CAS:1310-73-2	Protection from use of inorganic bases
P Hydrogen peroxide 30%	CAS: 7722-84-1	Protection from use of peroxides
T Formaldehyde 37%	CAS: 50-00-0	Protection from the use of aldehydes

RECOMENDATIONS OF USE Multi-purpose glove to protect the hand and/or the handled product, can also be used for general uses ■ The glove is suitable for contact with pure chemical products for a period of less than the permeation time indicated in the chemical resistance table ■ Results obtained under laboratory conditions may not reflect actual working time

WARNINGS This informative does not reflect real chemical resistance and penetration resistance which have been laboratory tested from samples taken from the palm only (except in cases where the glove is 400mm or more in length - where the cuff is also tested) and refer only to the tested sample. It may be different if the chemical is used in a mixture ■ It is recommended to check that the gloves are suitable for the intended use as the conditions at the workplace may differ from the test depending on the temperature, abrasion and degradation ■ If used, protective gloves may provide less resistance to the hazardous chemical due to changes in physical properties. The movements, entanglement, rubbing, deterioration caused by chemical contact, etc. can significantly reduce the actual usage time. For corrosive chemicals, degradation is the most important factor to consider when selecting chemical resistant gloves ■ Before using the product, try it in the real condition of use ■ The glove is a device subject to wear, so frequent replacement is recommended and in all cases where imperfections are observed ■ Wear gloves with dry, clean hands ■ In rare cases where transient hypersensitive skin reactions may occur, discontinue use and consult a doctor ■ The user must inform the manufacturer and the relative authority of any incidents occurring during the use of the protective equipment

STORAGE AND DISPOSAL Store the gloves in their original packaging in a dry place ■ Store at room temperature ■ Do not expose to direct sunlight and heat sources. ■ The expiry date is 3 years from the date of manufacture indicated on the sales package, and is valid for the product properly stored ■ Gloves must be disposed of according to local regulations ■ Declaration of Conformity on www.reflexx.com