

**i** **INFORMATIVE NOTE**

<b>Name: REFLEXX Care L550</b>	
<b>Date: 03/11/2017</b>	
<b>Rev: 00</b>	

**DESCRIPTION:** Disposable glove non-sterile  
**Material:** natural latex (protein content < 50 µg/g)  
**Shape:** ambidextrous; five fingers, beaded cuff  
**Color:** light yellow  
**Minimum length:** 240 mm  
**Thickness in the palm:** 0.10mm (+/- 0.03)  
**Finishing:** externally textured surface; internally powder free, chlorinated and treated with colloidal oats  
**Article/Sizes (CE EN 420):** L550/XS (6) – L550/S (7) – L550/M (8) – L550/L (9) – L550/XL (10)  
**Packaging:** 100 pieces / box; 10 boxes / carton  
**Tensile strength (455-2: 2009 / A1):** ≥ 6 N



**CONFORMITY AND PERFORMANCES**

**CE** Medical Device Class I Dir. 93/42/CEE – Dir 2007/47/CEE – EN 455 1,2,3, & 4  
 Personal Protective Equipment (PPE) cat I only for minor risks – Dir. 89/686/EEC - EN420 – EN374

**ISO 2859:** AQL 1,5 (G1 level)

**RECOMMENDATIONS OF USE**

- Multi-purpose glove to protect the hand and/or the product, it can be used for general purposes
- Powder-free glove, also recommended for people allergic and/or sensitive to corn starch
- Glove treated with colloidal oats, hydrating effect

**INSTRUCTIONS FOR USE**

- Before using the product, test it in real conditions of use
- Avoid contact with aggressive substances such as fats, oils, acids, hydrocarbons and other chemical agents
- The product contains latex rubber, that may cause allergic reactions of hypersensitivity (type I and IV)
- The glove is a device subject to wear and tear so it is recommended the frequent replacement , and in all cases where visually you can see imperfections.
- Wear gloves with dry clean hands
- Carefully read the Information Note before using the glove

**STORAGE AND DISPOSAL**

- Store the gloves in the original packaging in a dry place
- Avoid exposure to sunlight and heat sources
- Dispose the gloves in accordance with local regulations

