



## INFORMATIVE NOTE

**Name: FEATHER PLUS**

**Date: 18/10/2021 Rev: 09**

### MANUFACTURER

REFLEXX SpA  
Via Passeri, 2 46019 VIADANA (MN)  
Tel.: + 39 0375 758891 Fax + 39 0375 464504  
[info@reflexx.com](mailto:info@reflexx.com)  
[www.reflexx.com](http://www.reflexx.com)



### Imballo - Pack

100		10	
Taglia Size	Articolo Article	Taglia Size	Cartone Carton
FP S	8008990232157	S	8008990087013
FP M	8008990232144	M	8008990087006
FP L	8008990232151	L	8008990086993
FP XL	8008990232380	XL	8008990086058

**DESCRIPTION:** ■ Disposable powdered glove, not sterile ■ **Material:** natural rubber latex ■ **Shape:** ambidextrous, five fingers, with beaded cuff ■ **Colour:** white ■ **Feature:** smooth surface ■ **Packing:** 100 pieces/box; 10 boxes/carton ■ **ISO 2859:** AQL 1.5 (level G1) freedom from holes; AQL 2.5 (level G1) major defects; AQL 4 (level G1) minor defects

### CONFORMITY AND PERFORMANCES



cat I

**Personal Protective Equipment cat I** for minimal risks only  
Reg EU 2016/425  
EN 420:2003 + A1 2009



Reg CE 1935/2004 – Reg UE 10/2011 - DM 21/03/19. To such use carefully read the declaration of suitability for food contact (DICA).

DIMENSIONS		
Cod./Size	Length min.	Width (±5 mm)
FP/S (7)	240	85
FP/ M (8)	240	95
FP/ L (9)	240	106
FP/ XL (10)	240	116

<b>Thickness</b> mm	Palm (+/- 0.03)	0.07
---------------------	-----------------	------

<b>EN 455-2:2009/A1</b>	Tensile strength	≥6 N
-------------------------	------------------	------

**RECOMENDATIONS OF USE** Multi-purpose glove to protect the hand and/or the handled product ■ Avoid contact with aggressive substances such as fats, oils, acids, hydrocarbons and other chemical agents  
**WARNINGS** Before using the product, try it in the real condition of use ■ The glove is a device subject to wear, so frequent replacement is recommended and in all cases where imperfections are observed ■ Wear gloves with dry, clean hands ■ The product contains latex rubber. It may cause allergic reactions of type I hypersensitivity and delayed hypersensitivity type IV. In case of reaction, discontinue use and consult a doctor  
**STORAGE AND DISPOSAL** Store the gloves in their original packaging in a dry place. ■ Do not expose to direct sunlight and heat sources. ■ The expiry date is 3 years from the date of manufacture indicated on the sales package, and is valid for the product properly stored ■ Gloves must be disposed of according to local regulations ■ Product Declaration of Conformity on [www.reflexx.com](http://www.reflexx.com)