

INFORMATIVE NOTE

Name: REFLEXX 36		Date: 01/01/2024 Rev: 019																													
MANUFACTURER REFLEXX SpA Via Passeri, 2 46019 VIADANA (MN) Tel.: + 39 0375 758891 Fax + 39 0375 464504 info@reflexx.com www.reflexx.com		 <table><tr><th colspan="4">Imballo - Pack</th></tr><tr><th> 100 </th><th></th><th></th><th> 10 </th></tr><tr><th>Taglia Size</th><th>Articolo Article</th><th>Taglia Size</th><th>Cartone Carton</th></tr><tr><td>R36 S</td><td> 8032891631430</td><td>S</td><td> 8032891636435</td></tr><tr><td>R36 M</td><td> 8032891631447</td><td>M</td><td> 8032891636442</td></tr><tr><td>R36 L</td><td> 8032891631454</td><td>L</td><td> 8032891636459</td></tr><tr><td>R36 XL</td><td> 8032891631461</td><td>XL</td><td> 8032891636466</td></tr></table>		Imballo - Pack				100			10	Taglia Size	Articolo Article	Taglia Size	Cartone Carton	R36 S	 8032891631430	S	 8032891636435	R36 M	 8032891631447	M	 8032891636442	R36 L	 8032891631454	L	 8032891636459	R36 XL	 8032891631461	XL	 8032891636466
Imballo - Pack																															
100			10																												
Taglia Size	Articolo Article	Taglia Size	Cartone Carton																												
R36 S	 8032891631430	S	 8032891636435																												
R36 M	 8032891631447	M	 8032891636442																												
R36 L	 8032891631454	L	 8032891636459																												
R36 XL	 8032891631461	XL	 8032891636466																												

DESCRIPTION: ■ Disposable examination glove not sterile powder free ■ **Material:** Vinyl (PVC) ■ **Shape:** ambidextrous; five fingers, beaded cuff ■ **Colour:** white color ■ **Feature:** smooth surface ■ **Packing:** 100 pieces/box; 10 boxes/carton ■ **ISO 2859:** AQL 1.5 (level G1) freedom from holes; AQL 1.5 (level G1) for major defects and minor defects

CONFORMITY AND PERFORMANCES

Personal Protective Equipment cat III

0465 ANCI_CIMAC, Via Aguzzafame, 60/b 27029 Vigevano (PV)

Reg UE 2016/425 - EN ISO 374-1:2016 TYPE C - EN ISO 374-5:2016 VIRUS (ISO 16604:2004). - EN 374-2:2014 - EN 420:2003 + A1 2009

Medical Device

Reg UE 2017/745 – EN 455 1,2,3 & 4

(CND T01020201 NOT SURGICAL VINYL GLOVES)

- Complies with ASTM-F-1671 (viral penetration test) and ISO 10993-10 regulations (skin irritation and sensitization)



Reg CE 1935/2004 – Reg UE 10/2011 - DM 21/03/73. To such use carefully read the declaration of suitability for food contact (DICA).

DIMENSIONS		
Cod./Size	Length min.* (mm)	Width mm (±5)
R36/ S (7)	240	85
R36/ M (8)	240	95
R36/ L (9)	240	106
R36/ XL (10)	240	116
* Gloves for special applications 5.1.3		

Thickness mm	Palm (+/- 0.02)	0.07
--------------	-----------------	------

EN 455-2:2009/A1	Tensile strength	≥3.6 N
------------------	------------------	--------

REG EU 2016/425

EN ISO 374-1:2016 + A1:2018/Type C



EN ISO 374-5: 2016



EN ISO 374-1:2016				
Chemical compound :		EN 16523-1:2015 TYPE C		EN 374-4:2013
	COD	Performance level	Breakthrough time	% degradation
Sodium hydroxide 40%	K	2	>30'	52%
Hydrogen peroxide 30%	P	2	>30'	17%

RISKS FROM WHICH THE DPI IS DESIGNED TO PROTECT:

Code chemical compound	Protection from	
K Sodium hydroxide 40%	CAS:1310-73-2	Protection from use of inorganic bases
P Hydrogen peroxide 30%	CAS: 7722-84-1	Protection from use of peroxides

RECOMENDATIONS OF USE Multi-purpose glove to protect the hand and/or the handled product, can also be used for general uses ■ The glove is suitable for contact with pure chemical products for a

REFLEXX S.p.A. UNIPERSONALE Via PASSERI, 2 - 46019 VIADANA (MN) | TEL. + 39 0375 758891

Cap.Soc. 1.200.000 € (i.v.) | PIVA 02085450209 | M5UXCR1 | R.E.A. 223166 | info@reflexx.com www.reflexx.com



period of less than the permeation time indicated in the chemical resistance table ■ Results obtained under laboratory conditions may not reflect actual working time **WARNINGS** Before using the product, try it in the real condition of use ■ The glove is a device subject to wear, so frequent replacement is recommended and in all cases where imperfections are observed ■ Wear gloves with dry, clean hands ■ The product does not contain natural rubber latex ■ The product does not contain nickel ■ In case of transient skin reactions and hypersensitivity reactions, discontinue use and consult a doctor **STORAGE AND DISPOSAL** Store the gloves in their original packaging in a dry place ■ Do not expose to direct sunlight and heat sources. ■ The expiry date is 5 years from the date of manufacture indicated on the sales package, and is valid for the product properly stored ■ Gloves must be disposed of according to local regulations ■ Product Declaration of Conformity on www.reflexx.com