

 **INFORMATIVE NOTE**

<b>Name: REFLEXX 36</b>		<b>Date: 06/12/2024 Rev: 021</b>																									
<b>MANUFACTURER</b> REFLEXX SpA Via Passeri, 2 46019 VIADANA (MN) Tel.: + 39 0375 758891 Fax + 39 0375 464504 <a href="mailto:info@reflexx.com">info@reflexx.com</a> <a href="http://www.reflexx.com">www.reflexx.com</a>		<b>Imballo - Pack</b>  <table border="1"> <thead> <tr> <th colspan="2">100</th> <th colspan="2">10</th> </tr> <tr> <th>Taglia Size</th> <th>Articolo Article</th> <th>Taglia Size</th> <th>Cartone Carton</th> </tr> </thead> <tbody> <tr> <td>R36 S</td> <td>8032891631430</td> <td>S</td> <td>8032891636435</td> </tr> <tr> <td>R36 M</td> <td>8032891631447</td> <td>M</td> <td>8032891636442</td> </tr> <tr> <td>R36 L</td> <td>8032891631454</td> <td>L</td> <td>8032891636459</td> </tr> <tr> <td>R36 XL</td> <td>8032891631461</td> <td>XL</td> <td>8032891636466</td> </tr> </tbody> </table>		100		10		Taglia Size	Articolo Article	Taglia Size	Cartone Carton	R36 S	8032891631430	S	8032891636435	R36 M	8032891631447	M	8032891636442	R36 L	8032891631454	L	8032891636459	R36 XL	8032891631461	XL	8032891636466
100		10																									
Taglia Size	Articolo Article	Taglia Size	Cartone Carton																								
R36 S	8032891631430	S	8032891636435																								
R36 M	8032891631447	M	8032891636442																								
R36 L	8032891631454	L	8032891636459																								
R36 XL	8032891631461	XL	8032891636466																								

**DESCRIPTION:** ■ Disposable examination glove not sterile powder free ■ **Material:** Vinyl (PVC) ■ **Shape:** ambidextrous; five fingers, beaded cuff ■ **Colour:** white color ■ **Feature:** smooth surface ■ **Packing:** 100 pieces/box; 10 boxes/carton ■ **ISO 2859:** AQL 1.5 (level G1) freedom from holes; AQL 1.5 (level G1) for major defects and minor defects

**CONFORMITY AND PERFORMANCES**

**Personal Protective Equipment cat III**

0465 ANCI\_CIMAC, Via Aguzzafame, 60/b 27029 Vigevano (PV)

Reg EU 2016/425 - EN ISO 374-1:2016 TYPE C - EN ISO 374-5:2016 VIRUS (ISO 16604:2004). - EN 374-2:2014 - - EN 420:2003 + A1 2009 - EN 21420:2020

**Medical Device**

Reg UE 2017/745 – EN 455 1,2,3 & 4

(CND T01020201 NOT SURGICAL VINYL GLOVES)

- Complies with ASTM-F-1671 (viral penetration test) and ISO 10993-10 regulations (skin irritation and sensitization)



Reg CE 1935/2004 – Reg UE 10/2011 - DM 21/03/73. To such use carefully read the declaration of suitability for food contact (DICA).

DIMENSIONS		
Cod./Size	Length min.* (mm)	Width mm (±5)
R36/ S (7)	240	85
R36/ M (8)	240	95
R36/ L (9)	240	106
R36/ XL (10)	240	116

\* Gloves for special applications 5.1.3

<b>Thickness</b> mm	Palm (+/- 0.02)	0.07
------------------------	-----------------	------

<b>EN 455-2:2009/A1</b>	Tensile strength	≥3.6 N
-------------------------	------------------	--------

**REG EU 2016/425**

EN ISO 374-1:2016 + A1:2018/Type C



**EN ISO 374-5:2016**



EN ISO 374-1:2016				
Chemical compound :	EN 16523-1:2015 TYPE C			EN 374-4:2013
	COD	Performance level	Breakthrough time	% degradation
Sodium hydroxide 40%	K	6	>480'	3,4%
Hydrogen peroxide 30%	P	2	>30'	17,3%

RISKS FROM WHICH THE DPI IS DESIGNED TO PROTECT:

Code chemical compound	CAS	Protection from
K Sodium hydroxide 40%	CAS:1310-73-2	Protection from use of inorganic bases
P Hydrogen peroxide 30%	CAS: 7722-84-1	Protection from use of peroxides

**RECOMENDATIONS OF USE** Multi-purpose glove to protect the hand and/or the handled product, can also be used for general uses ■ The glove is suitable for contact with pure chemical products for a period of less than the permeation time indicated in the chemical resistance table ■ Results obtained under laboratory conditions may not reflect actual working time

**WARNINGS** Before using the product, try it in the real condition of use ■ The glove is a device subject to wear, so frequent replacement is recommended and in all cases where imperfections are observed ■ Wear gloves with dry, clean hands ■ The product does not contain natural rubber latex ■ The product does not contain nickel ■ In case of transient skin reactions and hypersensitivity reactions, discontinue use and consult a doctor ■ The user must inform the manufacturer and the relative authority of any incidents occurring during the use of the protective equipment

**STORAGE AND DISPOSAL** Store the gloves in their original packaging in a dry place ■ Do not expose to direct sunlight and heat sources. ■ The expiry date is 5 years from the date of manufacture indicated on the sales package, and is valid for the product properly stored ■ Gloves must be disposed of according to local regulations ■ Product Declaration of Conformity on [www.reflexx.com](http://www.reflexx.com)