


INFORMATIVE NOTE
Name: REFLEXX care L600
Date: 15/05/2024 Rev: 03
MANUFACTURER

 REFLEXX SpA
 Via Passeri, 2 46019 VIADANA (MN)
 Tel.: + 39 0375 758891 Fax + 39 0375 464504
info@reflexx.com
www.reflexx.com

Pack

Cod./Size	Box	Carton
L600/ XS (6)	8032891632130	8032891637135
L600/ S (7)	8032891632147	8032891637142
L600/ M (8)	8032891632154	8032891637159
L600/ L (9)	8032891632161	8032891637166
L600/ XL (10)	8032891632178	8032891637173

DESCRIPTION: ■ disposable examination glove powder free, not sterile ■ **Material:** natural latex (protein content < 50 µg/g) ■ **Shape:** ambidextrous; five fingers, beaded cuff ■ **Color:** green ■ **Feature:** internally chlorinated and treated with aloe vera, moisturizing effect; externally textured surface ■ **Packaging:** 100 pieces / box; 10 boxes / carton ■ **ISO 2859:** AQL 1.5 (Level G1) freedom from holes; AQL 2.5 (level G1) major defects; AQL 4 (level G1) minor defects

CONFORMITY AND PERFORMANCES
Personal Protective Equipment cat III

0465 ANCI_CIMAC, Via Aguzzafame, 60/b 27029 Vigevano (PV)

Reg EU 2016/425 - EN ISO 374-1:2016 TYPE C - EN ISO 374-5:2016 - EN 374-2:2014 - EN 420:2003 + A1 2009

Medical Device

Reg EU 2017/745 - EN 455 1,2,3 & 4

(CND T010201 NON SURGICAL LATEX GLOVES)

DIMENSIONS		
Cod./Size	Length. min* (mm)	Width. (±5 mm)
L600/ XS (6)	240	75
L600/ S (7)	240	85
L600/ M (8)	240	95
L600/ L (9)	240	106
L600/ XL (10)	240	116

* Gloves for special applications 5.1.3

Thickness mm	Palm (+/- 0.03)	0.12
------------------------	-----------------	------

EN 455-2:2009/A1	Tensile strength	≥6 N
-------------------------	------------------	------

EN ISO 374-5:2016	Protection from microorganisms (bacteria and fungi)	Virus protection ISO 16604:2004
	pass	pass

EN ISO 374-1:2016				
Chemical compound:	COD	EN 16523-1:2015 TYPE C		EN 374-4:2013
		Performance level	Breakthrough time	% degradation
Sodium hydroxide 40%	K	6	>480'	1.0%


RISKS FROM WHICH THE DPI IS DESIGNED TO PROTECT:

COD		Protection
K Sodium hydroxide 40%	CAS:1310-73-2	Protection from use of inorganic bases

RECOMENDATIONS OF USE Multi-purpose glove to protect the hand and/or the handled product, can also be used for general uses ■ The glove is suitable for contact with pure chemical products for a period of less than the permeation time indicated in the chemical resistance table ■ Results obtained under

laboratory conditions may not reflect actual working time **WARNINGS** Before using the product, try it in the real condition of use ■ The glove is a device subject to wear, so frequent replacement is recommended and in all cases where imperfections are observed ■ Wear gloves with dry, clean hands ■ The product contains natural latex which can cause allergic reactions of type I and IV hypersensitivity. In case of reaction, discontinue use and consult a doctor ■ The user must inform the manufacturer and the relative authority of any incidents occurring during the use of the protective equipment **STORAGE AND DISPOSAL** Store the gloves in their original packaging in a dry place. ■ Store at room temperature ■ Do not expose to direct sunlight and heat sources. ■ The expiry date is 3 years from the date of manufacture indicated on the sales package, and is valid for the product properly stored ■ Gloves must be disposed of according to local regulations ■ Declaration of Conformity on www.reflexx.com